



www.accutechortho.com

Richmond Location
12617 Wilfong Drive
Midlothian, VA 23112
Local: 804-378-3728
Toll Free: 1-800-413-4522
Fax: 804-794-2087

Suffolk Location
3345 Bridge Road, Suite 904
Suffolk, VA 23435
Local: 757-488-1919
Toll Free: 1-800-734-7855
Fax: 757-488-7340

ALL SHADED AREAS & RETURN DATE MUST BE FILLED OUT BEFORE WE CAN PROCESS YOUR ORDER

Patient's Name _____ Date Sent _____
(First) / (Last)
Drs. Name _____ Practice Name _____
Address _____ Office _____
City _____ State _____ Zip _____
Phone # _____
Doctor's Signature _____ License # _____

Return Date:

By 5:00 pm on this day

**Please keep
yellow copy for
your records**

Order Lab Supplies Here

Please put how many you would like to the right of each request

- ☐ Rxs _____
☐ Shipping Boxes _____
☐ Color Chart _____
☐ 2nd Day Shipping Labels _____
☐ Safety Bags _____
☐ Local white pickup Bags _____

- U L ClearBow Retainers**
☐ Circumferential (Wraparound)
☐ SLB (Soldered Labial Bow)

- U L Retainers**
☐ Hawley w/ Ball Clasps
☐ Hawley w/ "C" Clasps
☐ Hawley w/ Adams Clasp
☐ Hawley w/ Arrow Clasps
☐ Hawley w/ soldered "C" Clasps on 4's
☐ Hawley w/o Clasps
☐ Circumferential (Wraparound)
☐ SLB (Soldered Labial Bow)
☐ 3x3 Clip
☐ ACR-(Accutech Comfort Retainer)
☐ Anterior Bite Plate w/ Ball Clasps
☐ Anterior Bite Plate w/ Buccal Tube Clasps
☐ Essix (Invisible retainer)
☐ Tremont Cantilever Wrap Around Appl.

- Accessories**
☐ Labial acrylic
☐ Flat labial bow wire
☐ "z" spring
☐ "z" spring w/ Helix
☐ Add anterior bite plate
☐ Add plastic tooth/teeth- shade _____

- Spring Retainers**
☐ Modified Spring Retainer (good)
☐ Modified Spring Retainer +2 (better)
☐ Super + 2 Spring Retainer (best)
☐ Clasps Type _____
☐ Use TMA spring wire for best results
☐ Reset tooth/teeth # _____

- Distalizers**
☐ Distal Jet Molar Distalizer (Unilateral)
☐ Distal Jet Molar Distalizer (Bilateral)
☐ Pendulum ☐ Pend X
☐ T-Rex ☐ Penguin

- Fixed Appliances**
☐ 6 x 6 Lingual Arch with Loops
☐ 6 x 6 Lingual Arch without Loops
☐ Nance Appliance
☐ TPA
☐ Bonded 3 x 3 w/ tray
☐ Space Maintainer
☐ Habit Appliance
☐ Vertical - ☐ Horizontal
☐ BLR (Bonded Lingual Retainer)
Wire Type: ☐ Flat braided ☐ .018 Round
☐ Bond-Tic (Fixed anterior tooth replacement)
Tooth shade _____

- Herbts & Accessories**
☐ Accutech Standard Cantilever Design
☐ Design Your Own (Use special instructions)

- ☐ Use Hanks Hardware
☐ Crowns or Bands - You Provide & We Fit
☐ Crowns or Bands - We Provide & Fit
☐ Archwire Tubes
☐ Advancement Shims (1-5mm)
☐ R.P.E. Expansion Screw
☐ Soldered Wire Molar Rest ☐ Buttons
☐ Hooks ☐ Lingual Arch ☐ TPA

- Functionals**
☐ Basic Activator - Type _____
☐ Bionator 1 or 2
☐ Frankel - Type _____
☐ Twin Block
☐ Other _____

- Fixed Expanders**
☐ Standard R.P.E.
☐ Standard R.P.E. w/ Acrylic Pads
☐ Compact 2-Pin Mini-Screw R.P.E.
☐ SuperScrew™
☐ Haas R.P.E.
☐ Exspider
☐ Quad Helix
☐ Standard Bonded R.P.E.
☐ DeLuke Contoured Expander (DCE) -Bonded
☐ Freedom-Lock (Type) _____

**Plaster Study Models (From your impressions)
(6 Business Days in Lab)**

- ☐ Finished - *Angles, sculptured, soaped & polished*
☐ Unfinished - *Angles & sculptured*
☐ Rough cut - *Angles only*
☐ ABO models

- U L Splints & Options**
☐ Goldilocks ☐ Anterior Guidance
☐ Duraflex (Hard) ☐ Flat Plane
☐ Duraflex (Hard/Soft) ☐ No clasps
☐ All Hard Acrylic Splint ☐ Add ball clasps as indicated
(For all hard splints only)
☐ Use enclosed VO bite & do not change opening

- (Fabricated on the upper arch)*
☐ "Pro-Form" Athletic Mouth Guard w/ shock transmitter
☐ Standard Athletic Mouth Guard
☐ Add helmet strap

Mouth Guards available in colors listed below, ONLY

- ☐ Clear ☐ Blue ☐ Green ☐ Purple ☐ Yellow
☐ Orange ☐ Black ☐ White ☐ Red


SPECIAL INSTRUCTIONS

☐ **Please call me about this case**

Personalize Retainers with Acrylic colors, patterns &/or inlays

Max Color _____ Mand Color _____

Inlay picture (max Only) _____



All appliances returned to your office in 7 business days! GUARANTEED, or they're FREE!